

FERTILE®

MICROFLUIDIC SPERM SORTING CHIPS

Material Safety Data Sheet (MSDS)

Date updated August 20, 2016

Section 1 – Product and Company Identification

Product Name : FERTILE Microfluidic Sperm Sorting Chips
Company : Koek Biotechnology, Bioengineering and Medical Services Industry & Trade Ltd. Co.
Address : Dokuz Eylul University Inciralti Campus DEPART Health Olive Building Mithatpasa Street
No: 56/20-215 Balcova İzmir / TURKEY
Tel/Fax No. : +90 232 236 78 14
Charged by : Quality Assurance Department
Emergency contact No.: Same as above

Section 2 – Composition / Information on Ingredient

Single Product or Mixture : Mixture
Chemical Name : Poly(methyl methacrylate) (PMMA), Glass
CAS No. : PMMA 9011-14-7, Glass 1344-09-8

Section 3 – Dimension

Length: 76,40 mm
Width: 24,50 mm
Depth: 1,76 mm
Diameter of Outlets: 3,50 mm
Capacity of each Chamber: 15 µl

Section 4 – Sterilization

Gamma sterilized

Section 5 – Expiration Date

Expiration Date: 1 years after sterilization

Section 6 – Transport Information

It is not cold-storage products while delivery
Ensure no leak from container
Load with no harm, damage or falling.
Ensure no package collapsing.
Avoid direct sunlight.
It is not dangerous materials.
UN No. N/A